

Prevent spread of MDR pathogens...

With

SaniHigene

Surface, Aerial & Linen Disinfectant



Composition

Polymeric Biguanide Hydrochloride < 10%
 Alkyl Dimethyl Benzyl Ammonium Chloride & Didecyl Dimethyl Ammonium Chloride < 10%

Benefits

- Better than 5th generation QAC because of added PHMB
- Suitable for all surfaces, furniture, appliances, dental & general practice environment
- Aldehyde free, Corrosion free & Chlorine free formula
- Biodegradable

Microbial Efficacy

Bacteria

Mycobacterium tuberculosis

Fungi

Viruses

MDR pathogens: *Acinetobacter*, *Klebsiella spp.*, *Pseudomonas* & *Candida albicans*

Applications

EN 14476
 EN 13697
 prEN 16616
 Biodegradability Test

Area of Application	Method of Application	Dilution	Contact Time
Surface disinfection Critical areas	Mopping/wiping floors, walls, furniture, equipments, operation tables and OT lamps, etc.	0.5% (5 ml made up to 1000 ml with water) 2% (20 ml made up to 1000 ml with water) For heavy load	15 Min
Surface disinfection Non-critical areas	Mopping of waiting halls, corridors, reception areas, tables, chairs, table tops, wards, etc.	0.5% (5 ml made up to 1000 ml with water) 1% (10 ml made up to 1000 ml with water) For heavy load	5 Min-15 Min
Spill Management blood/ vomit /other body fluid	Pour solution on spills	10% (100 ml made up to 1000 ml with water)	5 Min
Fogging	To be used for every 1000 cu.ft. in high risk areas like OTs, ICUs, ICCUs, NICUs, Cath labs, AKDs, etc.	0.5% -1% (5ml -10 ml made up to 1000 ml with water)	20 Min fogging & 40 min OT closing

Parameter for Linen	Concentration	Soak Time	Temperature
Laundry / Washing Manual Cleaning	Non infected 0.5%	15 min	Room Temperature
	Infected 1%		
Automatic Cleaning	Non Infected 0.5%	10 min	70°C
	Infected 1%		

Packing
 1 L bottle and
 5 L jerry can



109112A4SCENGBAL00



SIRMAXO
 Chemicals Pvt. Ltd.
 Partner in infection control