

IRRIGATING CYSTOTOMES



D 5016	25G - 0.50 X 16 MM
D 5017	27G - 0.40 X 16 MM

Product Information :

Tip designed for opening of anterior capsule prior to phacoemulsification with a luer lock to allow irrigation for each cystome.
Ready for use in cataract surgery.

Instructions for use :

This anesthesia needle is designed for administration of topical anesthetics in ophthalmic surgery.

Packing :

Packed in a box of 10, each individually sterile in a pouch.

TARUN ENTERPRISES

H.O.: 8/8, Strachy Road, Prayagraj – 211001, U.P., INDIA
Phone: +91 - 8176080204, Email: info@optitecheyecare.com

B.O.: A96/1, Second Floor, Phase II, Mayapuri Industrial Area, New Delhi – 110064
Phone: 011-47070362, Email: delhi@optitecheyecare.com