

TEFLON BLOCK



The cornea is placed over the Teflon block and the donor button is then punched out.

Aspheric curvature to match corneal curvature..

Can be sterilized by dry heat, autoclave & E.T.O gas.

TARUN ENTERPRISES

H.O.: 8/8, Strachy Road, Prayagraj (Allahabad) – 211001, U.P., INDIA
Phone: +91 - 8176080204, Email: info@optitecheyecare.com

B.O.: A96/1, Second Floor, Phase II, Mayapuri Industrial Area, New Delhi – 110064
Phone: 011-47070362, Email: delhi@optitecheyecare.com