

## **DISPOSABLE CAPSULE HOOK**



Capsule hooks are used for temporary support of the capsular bag during a cataract surgery when the zonules are weak or in case of subluxated lenses.

- Capsule hooks are made of Nylon with an adjustable silicon lock
- Designed such that the hook is angled to be in the plane of the anterior capsule
- Easy to get a firm and reliable grip on the rim
- Tip of the hook is rounded to prevent any accidental damage to the peripheral capsule while placing the hooks in the bag
- Round tip also provides a larger support area per hook there by reducing the number of hooks required in each surgery
- Silicon lock helps in retaining the hook in retracted position during the surgery

**Packing :** Box pack of 5 Pcs., Sterile

## TARUN ENTERPRISES

**H.O.:** 8/8, Strachy Road, Prayagraj (Allahabad) – 211001, U.P., INDIA Phone: +91 - 8176080204, Email: info@optitecheyecare.com

🕮 www.optitecheyecare.com 🕓 +91 - 9335154556 🛛 🚹 OptitechOfficial

**B.O.:** A96/1, Second Floor, Phase II, Mayapuri Industrial Area, New Delhi – 110064 Phone: 011-47070362, Email: delhi@optitecheyecare.com