

CURVED HYDRODISSECTION CANNULA



D 5026

27G - 0.40 X 22 MM

Product Information:

The tip is designed to separate the cortex from the capsular bag. The flat end allows gentle insertion and suitable irrigation so as to separate the nucleus and the cortex from the capsule prior to phacoemulsification.

Ready for use in cataract surgery.

Instruction for use:

Hydrodissection cannulas are designed to dissect the crystalline lens nucleus.

Packaging :

Packed in a box of 50, each individually sterile in a pouch

TARUN ENTERPRISES

H.O.: 8/8, Strachy Road, Prayagraj – 211001, U.P., INDIA
Phone: +91 - 8176080204, Email: info@optitecheyecare.com

B.O.: A96/1, Second Floor, Phase II, Mayapuri Industrial Area, New Delhi – 110064
Phone: 011-47070362, Email: delhi@optitecheyecare.com



www.optitecheyecare.com



+91 - 9335154556



OptitechOfficial



[optitecheyecare](https://www.instagram.com/optitecheyecare)



[optitecheyecare](https://www.linkedin.com/company/optitecheyecare)