

## PERIBULBAR ANAESTHESIA NEEDLE



D 5001	23G - 0.60 X 22 MM
D 5002	25G - 0.50 X 22 MM
D 5003	27G - 0.40 X 22 MM

### Product Information :

Peribulbar, retrobulbar and peri-retrobulbar needles, ready for use in ophthalmic surgery.

### Instructions for use :

This anesthesia needle is designed for administration of topical anesthetics in ophthalmic surgery.

### Packing :

Packed in a box of 10, each individually sterile in a pouch.

---

### TARUN ENTERPRISES

H.O.: 8/8, Strachy Road, Prayagraj – 211001, U.P., INDIA  
Phone: +91 - 8176080204, Email: info@optitecheyecare.com

B.O.: A96/1, Second Floor, Phase II, Mayapuri Industrial Area, New Delhi – 110064  
Phone: 011-47070362, Email: delhi@optitecheyecare.com