

## SUCTION TREPHINE



A very useful instrument to cut button from the recipient cornea. The outer shell of the trephine holds the recipient cornea by vacuum and the internal blade is rotated to take out the diseased button.

Available in ready to use sterile packs.

### TARUN ENTERPRISES

**H.O.:** 8/8, Strachy Road, Prayagraj (Allahabad) – 211001, U.P., INDIA  
Phone: +91 - 8176080204, Email: [info@optitecheyecare.com](mailto:info@optitecheyecare.com)

**B.O.:** A96/1, Second Floor, Phase II, Mayapuri Industrial Area, New Delhi – 110064  
Phone: 011-47070362, Email: [delhi@optitecheyecare.com](mailto:delhi@optitecheyecare.com)